

Custom Solution - 27 components in HNO3 5%

Art. ID	LM24-CUS-70831
Unit	1 x 125 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 12 months (originally closed bottle), 6 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	µg/mL		
	Boron (B)	[7440-42-8]	10	µg/mL		
	Barium (Ba)	[7440-39-3]	10	µg/mL		
	Beryllium (Be)	[7440-41-7]	10	µg/mL		
	Cadmium (Cd)	[7440-43-9]	2	µg/mL		
	Cobalt (Co)	[7440-48-4]	5	µg/mL		
	Chromium (Cr)	[7440-47-3]	5	µg/mL		
	Copper (Cu)	[7440-50-8]	10	µg/mL		
	Iron (Fe)	[7439-89-6]	10	µg/mL		
	Lithium (Li)	[7439-93-2]	10	µg/mL		
	Manganese (Mn)	[7439-96-5]	5	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	10	µg/mL		
	Nickel (Ni)	[7440-02-0]	2	µg/mL		
	Lead (Pb)	[7439-92-1]	5	µg/mL		
	Tin (Sn)	[7440-31-5]	50	µg/mL		
	Strontium (Sr)	[7440-24-6]	10	µg/mL		
	Tellurium (Te)	[13494-80-9]	50	µg/mL		
	Titanium (Ti)	[7440-32-6]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	50	µg/mL		
	Vanadium (V)	[7440-62-2]	10	µg/mL		
	Zinc (Zn)	[7440-66-6]	10	µg/mL		
	Calcium (Ca)	[7440-70-2]	1500	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		
	Sodium (Na)	[7440-23-5]	1000	µg/mL		
	Potassium (K)	[7440-09-7]	200	µg/mL		
	Silicon (Si)	[7440-21-3]	250	µg/mL		
	Phosphorus (P)	[7723-14-0]	10	µg/mL		