

Copper base, low alloy, C Cu 2, Disc approx. 40 mm x 30 mm (check sample)

Art. ID SUS-C-Cu-2-30MM
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Please note that the values varies from lot to lot / Example values

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
approx. example value	Silver (Ag)	[7440-22-4]	0,5	%		
approx. example value	Nickel (Ni)	[7440-02-0]	1	%		
approx. example value	Gold (Au)	[7440-57-5]	0,005	%		
approx. example value	Lead (Pb)	[7439-92-1]	0,7	%		
approx. example value	Bismuth (Bi)	[7440-69-9]	0,04	%		
approx. example value	Antimony (Sb)	[7440-36-0]	0,6	%		
approx. example value	Chromium (Cr)	[7440-47-3]	<0,005	%		
approx. example value	Tin (Sn)	[7440-31-5]	0,4	%		
approx. example value	Iron (Fe)	[7439-89-6]	<0,01	%		
approx. example value	Zinc (Zn)	[7440-66-6]	0,1	%		
approx. example value	Manganese (Mn)	[7439-96-5]	<0,005	%		