

Copper base, low alloy, C Cu 3, Disc approx. 40 mm x 30 mm (check sample)

Art. ID SUS-C-Cu-3-30MM
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Please note that the values varies from lot to lot / Example values

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
approx. example value	Arsenic (As)	[7440-38-2]	0,0875	%		
approx. example value	Sulfur (S)	[7704-34-9]	0,0229	%		
approx. example value	Cadmium (Cd)	[7440-43-9]	0,0096	%		
approx. example value	Selenium (Se)	[7782-49-2]	0,0475	%		
approx. example value	Cobalt (Co)	[7440-48-4]	0,0496	%		
approx. example value	Tellurium (Te)	[13494-80-9]	0,0194	%		