

### **Cast Iron, truncated cone 34 x 15+mm**

Art. ID                      CKD-233  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,12	%		
	Manganese (Mn)	[7439-96-5]	0,26	%		
	Silicon (Si)	[7440-21-3]	2,59	%		
	Phosphorus (P)	[7723-14-0]	0,0033	%		
	Sulfur (S)	[7704-34-9]	0,14	%		
	Chromium (Cr)	[7440-47-3]	1,92	%		
	Nickel (Ni)	[7440-02-0]	0,062	%		
	Tin (Sn)	[7440-31-5]	0,033	%		
	Antimony (Sb)	[7440-36-0]	0,055	%		
	Bismuth (Bi)	[7440-69-9]	~0	%		
	Arsenic (As)	[7440-38-2]	0,073	%		
	Lead (Pb)	[7439-92-1]	0,027	%		
	Selenium (Se)	[7782-49-2]	0,018	%		
	Tellurium (Te)	[13494-80-9]	0	%		