

Cast Iron, truncated cone 35 x 15+mm

Art. ID CKD-234
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,46	%		
	Manganese (Mn)	[7439-96-5]	1,39	%		
	Silicon (Si)	[7440-21-3]	2,02	%		
	Phosphorus (P)	[7723-14-0]	0,38	%		
	Sulfur (S)	[7704-34-9]	0,068	%		
	Chromium (Cr)	[7440-47-3]	0,46	%		
	Nickel (Ni)	[7440-02-0]	0,305	%		
	Tin (Sn)	[7440-31-5]	0,026	%		
	Antimony (Sb)	[7440-36-0]	0,002	%		
	Bismuth (Bi)	[7440-69-9]	0	%		
	Arsenic (As)	[7440-38-2]	0,01	%		
	Lead (Pb)	[7439-92-1]	~0,002	%		
	Selenium (Se)	[7782-49-2]	0,006	%		
	Tellurium (Te)	[13494-80-9]	0,02	%		