

## **Cast Iron, truncated cone 35 x 15+mm**

Art. ID                      CKD-235  
Unit                            each  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,73	%		
	Manganese (Mn)	[7439-96-5]	1,86	%		
	Silicon (Si)	[7440-21-3]	0,92	%		
	Phosphorus (P)	[7723-14-0]	0,78	%		
	Sulfur (S)	[7704-34-9]	0,035	%		
	Chromium (Cr)	[7440-47-3]	0,41	%		
	Nickel (Ni)	[7440-02-0]	0,195	%		
	Tin (Sn)	[7440-31-5]	0,024	%		
	Antimony (Sb)	[7440-36-0]	0,016	%		
	Bismuth (Bi)	[7440-69-9]	0,008	%		
	Arsenic (As)	[7440-38-2]	0,1	%		
	Lead (Pb)	[7439-92-1]	0,011	%		
	Selenium (Se)	[7782-49-2]	0,0005	%		
	Tellurium (Te)	[13494-80-9]	0,033	%		