

Low Alloy Cast Iron, truncated pyramid 37 x 37 x 20mm

Art. ID CKD-243B
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,29	%		
	Manganese (Mn)	[7439-96-5]	0,466	%		
	Silicon (Si)	[7440-21-3]	2,44	%		
	Phosphorus (P)	[7723-14-0]	0,173	%		
	Sulfur (S)	[7704-34-9]	0,081	%		
	Nickel (Ni)	[7440-02-0]	0,098	%		
	Chromium (Cr)	[7440-47-3]	0,394	%		
	Copper (Cu)	[7440-50-8]	0,191	%		
	Molybdenum (Mo)	[7439-98-7]	0,252	%		
	Vanadium (V)	[7440-62-2]	0,227	%		
	Titanium (Ti)	[7440-32-6]	0,003	%		
	Aluminium (Al)	[7429-90-5]	0,013	%		
	Magnesium (Mg)	[7439-95-4]	0	%		
	Cerium (Ce)	[7440-45-1]	0	%		
	Boron (B)	[7440-42-8]	0,01	%		
	Tin (Sn)	[7440-31-5]	0,11	%		
	Antimony (Sb)	[7440-36-0]	0,079	%		
	Arsenic (As)	[7440-38-2]	0,078	%		
	Lead (Pb)	[7439-92-1]	0,013	%		
	Bismuth (Bi)	[7440-69-9]	0	%		
	Zinc (Zn)	[7440-66-6]	0,025	%		
	Selenium (Se)	[7782-49-2]	~0,016	%		
	Tellurium (Te)	[13494-80-9]	~0	%		
	Cobalt (Co)	[7440-48-4]	0,027	%		
	Tungsten (W)	[7440-33-7]	0,026	%		
	Niobium (Nb)	[7440-03-1]	0,024	%		
	Zirconium (Zr)	[7440-67-7]	0	%		
	Lanthanum (La)	[7439-91-0]	0	%		
	Nitrogen (N)	[7727-37-9]	0,005	%		
	Iron (Fe)	[7439-89-6]	~93	%		