

Custom Solution - 13 components in HNO3 5%, HF tr.

Art. ID	LM24-CUS-70860
Unit	1 x 125 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 8 months (originally closed bottle), 4 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Sulfur (S)	[7704-34-9]	250	mg/L		
	Boron (B)	[7440-42-8]	5	mg/L		
	Iron (Fe)	[7439-89-6]	5	mg/L		
	Manganese (Mn)	[7439-96-5]	5	mg/L		
	Silicon (Si)	[7440-21-3]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	0,5	mg/L		
	Copper (Cu)	[7440-50-8]	0,5	mg/L		
	Zinc (Zn)	[7440-66-6]	0,5	mg/L		