

Tool Steel / High-Speed / Certified Reference Material / AISI M-2 / T11302, Disc Ø 38 x 3 mm

Art. ID IARM-FeM2-18-X
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

The M2 grade, (IARM-FeM2-18) is a high speed tool steel. This W-Mo high speed steel grade has a relatively high C content, and good wear resistance.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Value	Aluminium (Al)	[7429-90-5]	0,014 ± 0,002	%		
Certified Value	Arsenic (As)	[7440-38-2]	0,008 ± 0,003	%		
Certified Value	Carbon (C)	[7440-44-0]	0,853 ± 0,006	%		
Certified Value	Cobalt (Co)	[7440-48-4]	0,28 ± 0,02	%		
Certified Value	Chromium (Cr)	[7440-47-3]	4,23 ± 0,07	%		
Certified Value	Copper (Cu)	[7440-50-8]	0,098 ± 0,003	%		
Certified Value	Iron (Fe)	[7439-89-6]	81,4 ± 0,5	%		
Certified Value	Manganese (Mn)	[7439-96-5]	0,337 ± 0,005	%		
Certified Value	Molybdenum (Mo)	[7439-98-7]	4,92 ± 0,05	%		
Certified Value	Nitrogen (N)	[7727-37-9]	0,0148 ± 0,0008	%		
Certified Value	Niobium (Nb)	[7440-03-1]	0,021 ± 0,005	%		
Certified Value	Nickel (Ni)	[7440-02-0]	0,182 ± 0,009	%		
Certified Value	Oxygen (O)	[7782-44-7]	0,0016 ± 0,0006	%		
Certified Value	Phosphorus (P)	[7723-14-0]	0,025 ± 0,002	%		
Certified Value	Sulfur (S)	[7704-34-9]	0,001 ± 0,0008	%		
Certified Value	Silicon (Si)	[7440-21-3]	0,26 ± 0,03	%		
Certified Value	Tin (Sn)	[7440-31-5]	0,007 ± 0,002	%		
Certified Value	Titanium (Ti)	[7440-32-6]	0,0016 ± 0,0007	%		
Certified Value	Vanadium (V)	[7440-62-2]	1,9 ± 0,02	%		
Certified Value	Tungsten (W)	[7440-33-7]	5,81 ± 0,06	%		
Indicative Value	Boron (B)	[7440-42-8]	~4	ppm		
Indicative Value	Bismuth (Bi)	[7440-69-9]	<50	ppm		
Indicative Value	Calcium (Ca)	[7440-70-2]	<50	ppm		
Indicative Value	Cadmium (Cd)	[7440-43-9]	<2	ppm		
Indicative Value	Hydrogen (H)	[1333-74-0]	<1	ppm		
Indicative Value	Magnesium (Mg)	[7439-95-4]	<10	ppm		
Indicative Value	Lead (Pb)	[7439-92-1]	~10	ppm		
Indicative Value	Antimony (Sb)	[7440-36-0]	~140	ppm		

Indicative Value	Selenium (Se)	[7782-49-2]	<50	ppm
Indicative Value	Tantalum (Ta)	[7440-25-7]	~23	ppm
Indicative Value	Zinc (Zn)	[7440-66-6]	~5	ppm
Indicative Value	Zirconium (Zr)	[7440-67-7]	~35	ppm