

Custom Solution - 11 components 1000mg/L each of in HNO3 5%, HF 0.5%

Art. ID LM24-CUS-70864
Unit 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Conservation 24 months. (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hafnium (Hf)	[7440-58-6]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Silicon (Si)	[7440-21-3]	1000	mg/L		
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Silver (Ag)	[7440-22-4]	1000	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Nickel (Ni)	[7440-02-0]	1000	mg/L		
	Magnesium (Mg)	[7439-95-4]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		