

Trace Elements in Water (lot 0523)

Art. ID NW-TMDA-62.3
Unit 250 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	117 ± 9	µg/L		
	Antimony (Sb)	[7440-36-0]	59,9 ± 4,3	µg/L		
	Arsenic (As)	[7440-38-2]	57,2 ± 4,7	µg/L		
	Barium (Ba)	[7440-39-3]	115 ± 7	µg/L		
	Beryllium (Be)	[7440-41-7]	54,3 ± 4,7	µg/L		
	Bismuth (Bi)	[7440-69-9]	49,0 ± 6,4	µg/L		
	Cadmium (Cd)	[7440-43-9]	91,7 ± 4,9	µg/L		
	Chromium (Cr)	[7440-47-3]	90,6 ± 5,5	µg/L		
	Cobalt (Co)	[7440-48-4]	92,0 ± 6,8	µg/L		
	Copper (Cu)	[7440-50-8]	87,0 ± 5,8	µg/L		
	Gallium (Ga)	[7440-55-3]	32,1 ± 1,4	µg/L		
	Iron (Fe)	[7439-89-6]	113 ± 8	µg/L		
	Lead (Pb)	[7439-92-1]	93,9 ± 6,4	µg/L		
	Lithium (Li)	[7439-93-2]	57,3 ± 5,7	µg/L		
	Manganese (Mn)	[7439-96-5]	91,4 ± 5,9	µg/L		
	Molybdenum (Mo)	[7439-98-7]	98,3 ± 6,2	µg/L		
	Nickel (Ni)	[7440-02-0]	92,4 ± 6,1	µg/L		
	Rubidium (Rb)	[7440-17-7]	15,6 ± 0,8	µg/L		
	Selenium (Se)	[7782-49-2]	52,9 ± 4,1	µg/L		
	Strontium (Sr)	[7440-24-6]	145 ± 10	µg/L		
	Thallium (Tl)	[7440-28-0]	50,8 ± 3,3	µg/L		
	Tin (Sn)	[7440-31-5]	104 ± 7	µg/L		
	Titanium (Ti)	[7440-32-6]	56,8 ± 3,5	µg/L		
	Uranium (U)	[7440-61-1]	54,0 ± 3,3	µg/L		
	Vanadium (V)	[7440-62-2]	112 ± 7	µg/L		
	Zinc (Zn)	[7440-66-6]	120 ± 10	µg/L		