

Low Alloy Steel, Certified Reference Material, Grade T23 / K40712, disc Ø 31 x 18 mm

Art. ID IARM-FET23-18-G
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	0,001	%		
	Carbon (C)	[7440-44-0]	0,068	%		
	Cobalt (Co)	[7440-48-4]	0,085	%		
	Chromium (Cr)	[7440-47-3]	2,47	%		
	Copper (Cu)	[7440-50-8]	0,046	%		
	Manganese (Mn)	[7439-96-5]	0,82	%		
	Molybdenum (Mo)	[7439-98-7]	0,261	%		
	Nitrogen (N)	[7727-37-9]	~0,003	%		
	Niobium (Nb)	[7440-03-1]	0,052	%		
	Nickel (Ni)	[7440-02-0]	0,53	%		
	Oxygen (O)	[7782-44-7]	0,0044	%		
	Phosphorus (P)	[7723-14-0]	0,012	%		
	Sulfur (S)	[7704-34-9]	0,006	%		
	Silicon (Si)	[7440-21-3]	0,18	%		
	Tin (Sn)	[7440-31-5]	~0,004	%		
	Titanium (Ti)	[7440-32-6]	~0,004	%		
	Vanadium (V)	[7440-62-2]	0,238	%		
	Tungsten (W)	[7440-33-7]	1,6	%		
	Boron (B)	[7440-42-8]	0,002	%		
	Carbon (C)	[7440-44-0]	0,068	%		
	Cobalt (Co)	[7440-48-4]	0,085	%		
	Chromium (Cr)	[7440-47-3]	2,47	%		
	Copper (Cu)	[7440-50-8]	0,046	%		
	Manganese (Mn)	[7439-96-5]	0,82	%		
	Molybdenum (Mo)	[7439-98-7]	0,261	%		
	Nitrogen (N)	[7727-37-9]	~0,003	%		
	Niobium (Nb)	[7440-03-1]	0,052	%		
	Nickel (Ni)	[7440-02-0]	0,53	%		
	Oxygen (O)	[7782-44-7]	0,0044	%		
	Phosphorus (P)	[7723-14-0]	0,012	%		
	Sulfur (S)	[7704-34-9]	0,006	%		

Silicon (Si)	[7440-21-3]	0,18	%
Tin (Sn)	[7440-31-5]	~0,004	%
Titanium (Ti)	[7440-32-6]	~0,004	%
Vanadium (V)	[7440-62-2]	0,238	%
Tungsten (W)	[7440-33-7]	1,6	%