

### **Custom Solution - 6 components 1000 mg/L each of in HNO3 2%**

Art. ID LM24-CUS-70895  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### **Description**

Conservation: 24 months (originally closed bottle), 12 months after opening.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Cerium (Ce)       | [7440-45-1] | 1000                | mg/L |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 1000                | mg/L |        |        |
|                  | Lithium (Li)      | [7439-93-2] | 1000                | mg/L |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 1000                | mg/L |        |        |
|                  | Thallium (Tl)     | [7440-28-0] | 1000                | mg/L |        |        |
|                  | Yttrium (Y)       | [7440-65-5] | 1000                | mg/L |        |        |