

Custom Solution - 12 components in HNO3 5%

Art. ID LM24-CUS-70900
Unit 1 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 12 months (originally closed bottle), 6 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Boron (B)	[7440-42-8]	5	µg/mL		
	Cadmium (Cd)	[7440-43-9]	5	µg/mL		
	Chromium (Cr)	[7440-47-3]	5	µg/mL		
	Copper (Cu)	[7440-50-8]	5	µg/mL		
	Nickel (Ni)	[7440-02-0]	5	µg/mL		
	Manganese (Mn)	[7439-96-5]	5	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	5	µg/mL		
	Zinc (Zn)	[7440-66-6]	5	µg/mL		
	Sodium (Na)	[7440-23-5]	500	µg/mL		
	Iron (Fe)	[7439-89-6]	500	µg/mL		
	Sulfur (S)	[7704-34-9]	500	µg/mL		
	Selenium (Se)	[7782-49-2]	2	µg/mL		