

Custom Solution - 21 components in Methanol

Art. ID LM24-CUS-70909
Unit 5 x 5 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Phenylphenol	[90-43-7]	10	µg/mL		
	2,3,4-Trichlorophenol	[15950-66-0]	10	µg/mL		
	2,3,5-Trichlorophenol	[933-78-8]	10	µg/mL		
	2,3,6-Trichlorophenol	[933-75-5]	10	µg/mL		
	2,4,5-Trichlorophenol	[95-95-4]	10	µg/mL		
	2,4,6-Trichlorophenol	[88-06-2]	10	µg/mL		
	3,4,5-Trichlorophenol	[609-19-8]	10	µg/mL		
	2,3,4,5-Tetrachlorophenol	[4901-51-3]	1	µg/mL		
	ol					
	2,3,4,6-Tetrachlorophenol	[58-90-2]	1	µg/mL		
	ol					
	2,3,5,6-Tetrachlorophenol	[935-95-5]	1	µg/mL		
	ol					
	Pentachlorophenol	[87-86-5]	1	µg/mL		
	2,3-Dichlorophenol	[576-24-9]	10	µg/mL		
	2,4-Dichlorophenol	[120-83-2]	10	µg/mL		
	2,5-Dichlorophenol	[583-78-8]	10	µg/mL		
	2,6-Dichlorophenol	[87-65-0]	10	µg/mL		
	3,4-Dichlorophenol	[95-77-2]	10	µg/mL		
	3,5-Dichlorophenol	[591-35-5]	10	µg/mL		
	2-Chlorophenol	[95-57-8]	10	µg/mL		
	3-Chlorophenol	[108-43-0]	10	µg/mL		
	4-Chlorophenol	[106-48-9]	10	µg/mL		
	Phenol	[108-95-2]	10	µg/mL		