

## **Custom Solution - 6 components in Water**

Art. ID LM24-CUS-70910  
Unit 50 x 500 mL bottle  
Deliverydetails No Dangerous Good /not restricted

### Description

Conservation: 8 months (originally closed bottle), 4 months after opening.

| Text/Information | Analyte/Parameter                          | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|------------|---------------------|------|--------|--------|
|                  | Nitrate (NO <sub>3</sub> -)                |            | 5                   | mg/L |        |        |
|                  | Nitrite (NO <sub>2</sub> -)                |            | 10                  | mg/L |        |        |
|                  | Chloride (Cl-)                             |            | 5                   | mg/L |        |        |
|                  | Sulfate (SO <sub>4</sub> <sup>2-</sup> )   |            | 25                  | mg/L |        |        |
|                  | Phosphate (PO <sub>4</sub> <sup>3-</sup> ) |            | 2                   | mg/L |        |        |
|                  | Fluoride (F-)                              |            | 0,8                 | mg/L |        |        |