

Custom Solution - 12 components 1000 µg/mL each of in Methanol

Art. ID LM24-CUS-70918
Unit 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phenol	[108-95-2]	1000	µg/mL		
	2-Chlorophenol	[95-57-8]	1000	µg/mL		
	3-Chlorophenol	[108-43-0]	1000	µg/mL		
	4-Chlorophenol	[106-48-9]	1000	µg/mL		
	4-Chloro-3-methylphenol	[59-50-7]	1000	µg/mL		
	2,4-Dichlorophenol	[120-83-2]	1000	µg/mL		
	2,4,5-Trichlorophenol	[95-95-4]	1000	µg/mL		
	Pentachlorophenol	[87-86-5]	1000	µg/mL		
	Bisphenol A	[80-05-7]	1000	µg/mL		
	4-tert-Octylphenol	[140-66-9]	1000	µg/mL		
	4-Nonylphenol (techn.)	[84852-15-3]	1000	µg/mL		
	2,4,6-Trichlorophenol	[88-06-2]	1000	µg/mL		