

Custom Solution - 9 components 100 mg/L each of in HNO3 0.2%, HCl 0.1%

Art. ID LM24-CUS-70919
Unit 1 x 250 mL bottle
Deliverydetails No Dangerous Good /not restricted

Description

Conservation: 8 months (originally closed bottle), 4 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Boron (B)	[7440-42-8]	1	mg/L		
	Silicon (Si)	[7440-21-3]	1	mg/L		
	Copper (Cu)	[7440-50-8]	1	mg/L		
	Zinc (Zn)	[7440-66-6]	0,5	mg/L		
	Magnesium (Mg)	[7439-95-4]	20	mg/L		
	Sodium (Na)	[7440-23-5]	25	mg/L		
	Potassium (K)	[7440-09-7]	5	mg/L		
	Iron (Fe)	[7439-89-6]	0,1	mg/L		