

## **Custom Solution - 4 components, 5 mg/L each in Toluene**

Art. ID LM24-CUS-70924  
Unit 5 x 1 mL  
Deliverydetails De Minimis UN1294 II E2

### Description

Conservation 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Phthalic acid,bis-2-ethylhexylester	Bis(2-ethylhexyl)phthalate	[117-81-7]	5	mg/L		
Phthalic acid, benzylbutyl ester	Benzylbutyl phthalate	[85-68-7]	5	mg/L		
Phthalic acid, bis-butyl ester	Di-n-butyl phthalate	[84-74-2]	5	mg/L		
Phthalic acid, bis-iso-butyl ester	Phthalic acid, bis-iso-butyl ester	[84-69-5]	5	mg/L		