

Custom Solution - 7 components 2000 mg/L in Ethanol

Art. ID LM24-CUS-70925
Unit 3 x 10 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1170 III E1

Description

Conservation: 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Safrole	[94-59-7]	2000	mg/L		
	Isosafrole	[120-58-1]	2000	mg/L		
	(-)-Fenchone	[7787-20-4]	2000	mg/L		
	4-Allylanisole	[140-67-0]	2000	mg/L		
	alpha-Thujone	[546-80-5]	2000	mg/L		
	Eugenol methyl ether	[93-15-2]	2000	mg/L		
	Pulegone	[89-82-7]	2000	mg/L		