

Custom solution - 25 components 100 mg/L each of in Acetone

Art. ID LM24-CUS-70926
Unit 6 x 1 mL (ampoule)
Deliverydetails De Minimis UN1090 II E2

Description

Conservation: 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alachlor	[15972-60-8]	100	mg/L		
	Ametryn	[834-12-8]	100	mg/L		
	Atrazine	[1912-24-9]	100	mg/L		
	Atrazine-desethyl	[6190-65-4]	100	mg/L		
	Terbutylazine-desethyl	[30125-63-4]	100	mg/L		
	Atrazine-desisopropyl	[1007-28-9]	100	mg/L		
	Desmetryn	[1014-69-3]	100	mg/L		
	Dimetridazole	[551-92-8]	100	mg/L		
	Bromacil	[314-40-9]	100	mg/L		
	Cyanazine	[21725-46-2]	100	mg/L		
	Metolachlor	[51218-45-2]	100	mg/L		
	Molinate	[2212-67-1]	100	mg/L		
	Oxadiazon	[19666-30-9]	100	mg/L		
	Prometryn	[7287-19-6]	100	mg/L		
	Propazine	[139-40-2]	100	mg/L		
	Sebutylazine	[7286-69-3]	100	mg/L		
	Simazine	[122-34-9]	100	mg/L		
	Terbutylazine	[5915-41-3]	100	mg/L		
	Terbutryn	[886-50-0]	100	mg/L		
	Tris(2-chloroethyl)phosphate	[115-96-8]	100	mg/L		
	2,6-Dichlorobenzamide	[2008-58-4]	100	mg/L		
	Tris(2-chloroisopropyl)phosphate	[13674-84-5]	100	mg/L		
	Pentetrazol	[54-95-5]	100	mg/L		
	Bentazone	[25057-89-0]	100	mg/L		
	Atrazine-desethyl-desisopropyl	[3397-62-4]	100	mg/L		