

## **Stainless Steel, Chips 150 g (JAPAN)**

Art. ID JSS-610-11  
Unit 150 g (chips)  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Carbon (C)        | [7440-44-0] | 1,252               | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 9,72                | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 4,04                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,0239              | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,301               | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 3,05                | %    |        |        |
|                  | Nitrogen (N)      | [7727-37-9] | 0,0349              | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,046               | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,0171              | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,00014             | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,31                | %    |        |        |
|                  | Vanadium (V)      | [7440-62-2] | 3,32                | %    |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 9,22                | %    |        |        |