

Stainless Steel, Chips 150 g **Last of stock 01/11 (JAPAN)

Art. ID JSS-670-3
Unit 150 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	~0,055	%		
	Carbon (C)	[7440-44-0]	0,008	%		
	Chromium (Cr)	[7440-47-3]	11,23	%		
	Copper (Cu)	[7440-50-8]	0,044	%		
	Manganese (Mn)	[7439-96-5]	0,53	%		
	Molybdenum (Mo)	[7439-98-7]	0,02	%		
	Nitrogen (N)	[7727-37-9]	0,0064	%		
	Nickel (Ni)	[7440-02-0]	0,21	%		
	Phosphorus (P)	[7723-14-0]	0,02	%		
	Sulfur (S)	[7704-34-9]	~0,0005	%		
	Silicon (Si)	[7440-21-3]	0,55	%		
	Titanium (Ti)	[7440-32-6]	0,26	%		