

Custom Solution - 9 components 1000 mg/L each of in 5% HNO3

Art. ID LM24-CUS-70933
Unit 125 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation: 24 months (originally closed bottle), 12 month after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Boron (B)	[7440-42-8]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Tellurium (Te)	[13494-80-9]	1000	mg/L		
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Selenium (Se)	[7782-49-2]	1000	mg/L		
	Vanadium (V)	[7440-62-2]	1000	mg/L		
	Beryllium (Be)	[7440-41-7]	1000	mg/L		
	Mercury (Hg)	[7439-97-6]	1000	mg/L		