

Custom Solution - Fluoranthene D10 100 mg/L in Cyclohexane

Art. ID LM24-CUS-70943
Unit 1 x 1 mL ampoule
Deliverydetails De Minimis UN1145 II E2

Description

Conservation: 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Fluoranthene D10 | [93951-69-0] | 100 | mg/L | | |