

### **Custom Solution - 13 components in HNO3 5%**

Art. ID LM24-CUS-70950  
Unit 125 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### **Description**

Conservation: 24 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	500	mg/L		
	Beryllium (Be)	[7440-41-7]	10	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Potassium (K)	[7440-09-7]	10000	mg/L		
	Strontium (Sr)	[7440-24-6]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Titanium (Ti)	[7440-32-6]	50	mg/L		
	Molybdenum (Mo)	[7439-98-7]	50	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Tin (Sn)	[7440-31-5]	500	mg/L		