

Custom Solution - 2 components in n-Hexane

Art. ID LM24-CUS-70977
Unit 2 x 1000 mL (bottle)
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1208 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	17500	mg/L		
	Styrene	[100-42-5]	2500	mg/L		