

Custom Solution - 37 components in Methanol

Art. ID LM24-CUS-70978
Unit 5 x 9 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl-tert.butylether (MTBE)	[1634-04-4]	100	mg/L		
	cis-1,2-Dichloroethene	[156-59-2]	100	mg/L		
	Chloroform (Trichlorome thane)	[67-66-3]	100	mg/L		
	1,1,1-Trichloroethane	[71-55-6]	100	mg/L		
	Carbon tetrachloride (T etrachloromethane)	[56-23-5]	100	mg/L		
	Trichloroethene	[79-01-6]	100	mg/L		
	Tetrachloroethene	[127-18-4]	100	mg/L		
	1,2-Dichloroethane	[107-06-2]	100	mg/L		
	Bromoform (Tribromometh ane)	[75-25-2]	100	mg/L		
	Dibromochloromethane	[124-48-1]	100	mg/L		
	Bromodichloromethane	[75-27-4]	100	mg/L		
	Benzene	[71-43-2]	100	mg/L		
	Toluene	[108-88-3]	100	mg/L		
	Ethylbenzene	[100-41-4]	100	mg/L		
	o-Xylene	[95-47-6]	100	mg/L		
	m-Xylene	[108-38-3]	50	mg/L		
	p-Xylene	[106-42-3]	50	mg/L		
	Naphthalene	[91-20-3]	100	mg/L		
	1,1-Dichloroethene	[75-35-4]	100	mg/L		
	1,1-Dichloroethane	[75-34-3]	100	mg/L		
	trans-1,2-Dichloroethen e	[156-60-5]	100	mg/L		
	Vinylchloride	[75-01-4]	100	mg/L		
	Styrene	[100-42-5]	100	mg/L		
	Methylene chloride (Dic	[75-09-2]	100	mg/L		

chloromethane)			
1,1,2-Trichloroethane	[79-00-5]	100	mg/L
1,1,2,2-Tetrachloroethane	[79-34-5]	100	mg/L
1,1,1,2-Tetrachloroethane	[630-20-6]	100	mg/L
3-Ethyltoluene	[620-14-4]	100	mg/L
1,3,5-Trimethylbenzene	[108-67-8]	100	mg/L
1,2,4-Trimethylbenzene	[95-63-6]	100	mg/L
Acrylic acid-methyl ester	[96-33-3]	1000	mg/L
Methacrylic acid-methyl ester	[80-62-6]	1000	mg/L
1,2-Dibromoethane	[106-93-4]	100	mg/L
1,4-Dioxane	[123-91-1]	1000	mg/L
tert.-Butanol	[75-65-0]	100	mg/L
Formic acid-tert-butyl ester	[762-75-4]	100	mg/L
5-Methyl-2-hexanone	[110-12-3]	100	mg/L