

### **Total Phenolics (4-AAP), WS QR**

Art. ID                      ERA-535QR  
Unit                         2 mL

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4-Dichlorophenol	[120-83-2]	0,0052-0,5	mg/L		
	2-Chlorophenol	[95-57-8]	0,0045-05	mg/L		
	Phenol	[108-95-2]	0,004-0,5	mg/L		
	Phenolics, total		0,001-0,5	mg/L		