

Total Phenolics (4-AAP), WS QR

Art. ID ERA-535QR
Unit 2 mL

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|------|--------|--------|
| | 2,4-Dichlorophenol | [120-83-2] | 0,0052-0,5 | mg/L | | |
| | 2-Chlorophenol | [95-57-8] | 0,0045-05 | mg/L | | |
| | Phenol | [108-95-2] | 0,004-0,5 | mg/L | | |
| | Phenolics, total | | 0,001-0,5 | mg/L | | |