

Ultra Low-Level TOC Standards Kit - Anatel 643

Art. ID ERA-14201
Unit kit

Description

Ultra Low-Level NIST Sucrose CRM -- This Kit contains 7 reference materials including 1 blank. The levels are 0.050, 0.060, 0.070, 0.080, 0.090, 0.100 mg/L. Standards are also available individually from Catalog (i.e., see p/n 13300, 13350, 13351, 13352, 13353, 13354, 13355). Unpreserved NIST Sucrose standards. Appropriate for Anatel 643 or A-1000. Do not μ Se for Sievers 400 or Thornton. 30-day Expiration