

PAT700 Calibration and Validation Control Kit

Art. ID ERA-19225
Unit Kit

Description

Set of 10 bottles includes (1) 100 µS/cm Conductivity prepared from KCl, (1) Reagent Water Rw, (1) 0.500 mg/L C µSP Sucrose Rs, (1) 0.500 mg/L C 1,4-Benzoquinone Rss, (1) Blank, (1) 0.250 mg/L C NIST Sucrose, (1) .500 mg/L C NIST Sucrose, (1) 0.750 mg/L C NIST Sucrose, (1) Blank, (1) 0.500 mg/L C NIST Sucrose /// Waters ERA will continue to guarantee a minimum shelf life of 30 days from date of shipment