

## **Sievers Model 800 Low Level Validation Set**

Art. ID                      ERA-19520  
Unit                         Set

### Description

Set includes 37 vials, (10) blanks, (4) 25.0 mg/L C NIST KHP, (4) 25.0 mg/L IC NIST NaHCO<sub>3</sub>, (2) 0.250 mg/L C NIST KHP, (9) 0.500 mg/L C NIST KHP, (2) 0.750 mg/L C NIST KHP, and (2) µSP System Suitability Sets.  
Replaces Sievers® Part STD 38101.