

TOC Calibration Kit for A-J N/C 2100 & 3100 -250mL

Art. ID ERA-19970
Unit Kit

Description

Set includes eight 250 mL amber bottles, (1) Water blank, (1) 1.00 mg/L C NIST KHP, (1) 2.50 mg/L C NIST KHP, (1) 5.00 mg/L C NIST KHP, (1) 10.0 mg/L C NIST KHP, (1) 25.0 mg/L C NIST KHP, (1) 50.0 mg/L C NIST KHP, (1) 100 mg/L C NIST KHP.