

PAT700 JP Calibration and Validation Control Kit

Art. ID ERA-19225J
Unit Kit

Description

Set of 9 bottles includes (1) 100 µS/cm Conductivity prepared from KCl, (3) Blanks, (1) 0.500 mg/L Sodium Dodecylbenzene Sulfonate (SDBS) Rss, (1) 0.500 mg/L C NIST KHP, (2) 0.750 mg/L C NIST KHP, (1) 1.250 mg/L C NIST KHP with RFID Tags for JP Compliance. µSe for PAT700