

Custom Solution - 3 components 2000 mg/L each of in Methanol

Art. ID LM24-CUS-70979
Unit 4 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acetone	[67-64-1]	2000	mg/L		
	n-Hexane	[110-54-3]	2000	mg/L		
	2-Propanol	[67-63-0]	2000	mg/L		