

### **Custom Solution - 11 components 100 mg/L each of in Methanol**

Art. ID LM24-CUS-70982  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,3,4-Tetrachlorobenzene	[634-66-2]	100	mg/L		
	1,2,3,5-Tetrachlorobenzene	[634-90-2]	100	mg/L		
	4-Chloro-3-methylphenol	[59-50-7]	100	mg/L		
	Ethyleneglycol-monomethyl ether	[109-86-4]	100	mg/L		
	Dimethyl-sulfoxide	[67-68-5]	100	mg/L		
	4-Methyl-2-pentanone	[108-10-1]	100	mg/L		
	2,4-Dinitrotoluene	[121-14-2]	100	mg/L		
	2-Propanol	[67-63-0]	100	mg/L		
	Propionitrile	[107-12-0]	100	mg/L		
	Ethyleneglycol-monoethyl ether	[110-80-5]	100	mg/L		
	1-Propanol	[71-23-8]	100	mg/L		