

Custom Solution - 4 components 1000 µg/mL each of in Methanol

Art. ID LM24-CUS-70983
Unit 1 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methyl-tert.butylether (MTBE)	[1634-04-4]	1000	µg/mL		
	Ethyl-tert-butylether (ETBE)	[637-92-3]	1000	µg/mL		
	tert-Amyl-methyl ether	[994-05-8]	1000	µg/mL		
	Diisopropylether	[108-20-3]	1000	µg/mL		