

Custom Solution - 30 components 1000 mg/L each of in Methanol

Art. ID LM24-CUS-70984
Unit 3 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ethyleneglycol-monoethyl ether	[110-80-5]	1000	mg/L		
	Ethyleneglycol-monobutyl ether	[111-76-2]	1000	mg/L		
	Ethyleneglycol-monomethyl ether	[109-86-4]	1000	mg/L		
	Ethanol	[64-17-5]	1000	mg/L		
	tert.-Butanol	[75-65-0]	1000	mg/L		
	1-Butanol	[71-36-3]	1000	mg/L		
	2-Butanol	[78-92-2]	1000	mg/L		
	2-Propanol	[67-63-0]	1000	mg/L		
	1-Propanol	[71-23-8]	1000	mg/L		
	Allyl alcohol	[107-18-6]	1000	mg/L		
	Acrylic acid-ethyl ester	[140-88-5]	1000	mg/L		
	Butylacetate	[123-86-4]	1000	mg/L		
	Acetic acid-methyl ester	[79-20-9]	1000	mg/L		
	Acrylic acid-methyl ester	[96-33-3]	1000	mg/L		
	Ethylacetate	[141-78-6]	1000	mg/L		
	Acetone	[67-64-1]	1000	mg/L		
	2-Butanone	[78-93-3]	1000	mg/L		
	2-Hexanone	[591-78-6]	1000	mg/L		
	4-Methyl-2-pentanone	[108-10-1]	1000	mg/L		
	2-Pentanone	[107-87-9]	1000	mg/L		
	Tetrahydrofuran	[109-99-9]	1000	mg/L		
	Ethylene oxide	[75-21-8]	1000	mg/L		
	1,3-Butadiene	[106-99-0]	1000	mg/L		

Cyclohexane	[110-82-7]	1000	mg/L
Acetonitrile	[75-05-8]	1000	mg/L
Sulfolan	[126-33-0]	1000	mg/L
1,4-Dioxane	[123-91-1]	1000	mg/L
Diethylether	[60-29-7]	1000	mg/L
Methyl-tert.butylether (MTBE)	[1634-04-4]	1000	mg/L
Ethyl-tert-butylether (ETBE)	[637-92-3]	1000	mg/L