

PCB 69 500 µg/mL in Iso-octane

Art. ID LM24-OP-4010801
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,3',4,6-Tetrachlorobiphenyl (PCB 69) | [60233-24-1] | 500 | µg/mL | | |