

PCB 166 500 µg/mL in Iso-octane

Art. ID LM24-OP-4010905
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|--------------|---------------------|-------|--------|--------|
| | 2,3,4,4',5,6-Hexachloro biphenyl (PCB 166) | [41411-63-6] | 500 | µg/mL | | |