

PCB 173 100 µg/mL in Iso-octane

Art. ID LM24-OP-4010910
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Validity: 24 months

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|--------------|---------------------|-------|--------|--------|
| | 2,2',3,3',4,5,6-Heptachlorobiphenyl (PCB 173) | [68194-16-1] | 100 | µg/mL | | |