

ASTM Method D2887/ ISO 3294 Calibration Mix (set of 5x1 mL) - 20 components 0,1 % (w/w) each of in Chloroform

Art. ID LM24-OP-4011604
Unit 5x1 mL
Deliverydetails De Minimis UN1888 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Decane	[124-18-5]	0,1	%		
	n-Dodecane	[112-40-3]	0,1	%		
	n-Dotriacontane	[544-85-4]	0,1	%		
	n-Eicosane	[112-95-8]	0,1	%		
	n-Heptadecane	[629-78-7]	0,1	%		
	n-Heptane	[142-82-5]	0,1	%		
	n-Hexadecane	[544-76-3]	0,1	%		
	n-Hexane	[110-54-3]	0,1	%		
	n-Hexatriacontane	[630-06-8]	0,1	%		
	n-Nonane	[111-84-2]	0,1	%		
	n-Octacosane	[630-02-4]	0,1	%		
	n-Octadecane	[593-45-3]	0,1	%		
	n-Octane	[111-65-9]	0,1	%		
	n-Pentadecane	[629-62-9]	0,1	%		
	n-Pentane	[109-66-0]	0,1	%		
	n-Tetracontane	[4181-95-7]	0,1	%		
	n-Tetracosane	[646-31-1]	0,1	%		
	n-Tetradecane	[629-59-4]	0,1	%		
	n-Tetratetracontane	[7098-22-8]	0,1	%		
	n-Undecane	[1120-21-4]	0,1	%		