

**ASTM Method D5580 Aromatics Quantitative Calibration Mixture with IS (set of 10x1 mL):
7 components**

Art. ID LM24-OP-4011609

Unit 10x1 mL

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	iso-Octane	[540-84-1]	72,3	%		
	Benzene	[71-43-2]	1,9	%		
	Toluene	[108-88-3]	2,3	%		
	2-Hexanone	[591-78-6]	10	%		
	Ethylbenzene	[100-41-4]	4,5	%		
	o-Xylene	[95-47-6]	4,5	%		
	1,2,4-Trimethylbenzene	[95-63-6]	4,5	%		