

ASTM Method D5580 Aromatics Quantitative Calibration Mixture with IS: 7 components

Art. ID LM24-OP-4011610

Unit 1 mL

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	iso-Octane	[540-84-1]	72,8	%		
	Benzene	[71-43-2]	4,5	%		
	Toluene	[108-88-3]	1	%		
	2-Hexanone	[591-78-6]	10	%		
	Ethylbenzene	[100-41-4]	9	%		
	o-Xylene	[95-47-6]	0,5	%		
	1,2,4-Trimethylbenzene	[95-63-6]	2,2	%		