

ASTM Method D7423 Oxygenates Calibration Mixture (set of 10x2 mL): 24 components in n-Hexane

Art. ID LM24-OP-4011611
Unit 10x2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1208 II E2

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1-Propanol	[71-23-8]	10	mg/L		
	2-Butanol	[78-92-2]	10	mg/L		
	2-Butanone	[78-93-3]	10	mg/L		
	Isovaleraldehyde	[590-86-3]	10	mg/L		
	Acetaldehyde	[75-07-0]	10	mg/L		
	Acetone	[67-64-1]	10	mg/L		
	Allyl alcohol	[107-18-6]	10	mg/L		
	Butyraldehyde	[123-72-8]	10	mg/L		
	1-Butanol	[71-36-3]	10	mg/L		
	Diethylether	[60-29-7]	10	mg/L		
	Diisopropylether	[108-20-3]	10	mg/L		
	Dimethyl Ether	[115-10-6]	10	mg/L		
	Ethanol	[64-17-5]	10	mg/L		
	Ethyl-tert-butylether (ETBE)	[637-92-3]	10	mg/L		
	2-Methyl-1-propanol	[78-83-1]	10	mg/L		
	Isobutyraldehyde	[78-84-2]	10	mg/L		
	2-Propanol	[67-63-0]	10	mg/L		
	Methanol	[67-56-1]	10	mg/L		
	Methyl-tert.butylether (MTBE)	[1634-04-4]	10	mg/L		
	Valeraldehyde	[110-62-3]	10	mg/L		
	Propionaldehyde	[123-38-6]	10	mg/L		
	Propyl Ether	[111-43-3]	10	mg/L		
	tert-Amyl-methyl ether	[994-05-8]	10	mg/L		
	tert.-Butanol	[75-65-0]	10	mg/L		