

## **ASTM Method D7845 Calibration Level 4: 34 components in Toluene**

Art. ID LM24-OP-4011630  
Unit 1 mL  
Deliverydetails De Minimis UN1294 II E2

### Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(R)-(+)-Limonene	[5989-27-5]	150	mg/kg		
	(+)-beta-Pinene	[19902-08-0]	150	mg/kg		
	1-Methylnaphthalene	[90-12-0]	240	mg/kg		
	1-Phenoxy-2-propanol	[770-35-4]	500	mg/kg		
	1-Phenylethanol	[98-85-1]	240	mg/kg		
	2,4-Dimethylphenol	[105-67-9]	240	mg/kg		
	2,4-Dimethylstyrene	[2234-20-0]	240	mg/kg		
	2,5-Dimethylstyrene	[2039-89-6]	240	mg/kg		
	2-Ethylphenol	[90-00-6]	500	mg/kg		
	2-Methylnaphthalene	[91-57-6]	240	mg/kg		
	2-Methylstyrene	[611-15-4]	150	mg/kg		
	2-Phenoxypropanol	[4169-04-4]	500	mg/kg		
	2-Phenylethanol	[60-12-8]	240	mg/kg		
	3-Ethylphenol	[620-17-7]	240	mg/kg		
	3-Methylstyrene	[100-80-1]	150	mg/kg		
	4-Ethylphenol	[123-07-9]	240	mg/kg		
	4-Isopropylphenol	[99-89-8]	500	mg/kg		
	4-Methylstyrene	[622-97-9]	150	mg/kg		
	Acrylic acid-butyl ester	[141-32-2]	150	mg/kg		
	r					
	alpha-Methylstyrene	[98-83-9]	150	mg/kg		
	1-Butanol	[71-36-3]	150	mg/kg		
	Cyclohexanol	[108-93-0]	150	mg/kg		
	Di-n-butylether	[142-96-1]	150	mg/kg		
	Dicyclopentadiene	[77-73-6]	150	mg/kg		
	alpha-Pinene	[80-56-8]	150	mg/kg		
	Ethylbenzene-D10	[25837-05-2]	500	mg/kg		
	Indene	[95-13-6]	240	mg/kg		
	Naphthalene	[91-20-3]	240	mg/kg		

p,alpha-Dimethylstyrene	[1195-32-0]	240	mg/kg
Phenol	[108-95-2]	240	mg/kg
Ethyleneglycol-monophen yl ether	[122-99-6]	500	mg/kg
Styrene	[100-42-5]	150	mg/kg
1-Phenyl-1,2-ethanediol	[93-56-1]	240	mg/kg
trans-beta-Methylstyren e	[873-66-5]	150	mg/kg