

ASTM Method D2887 Hydrocarbon Window Defining Solution: 24 components 500 mg/L each of in Chloroform

Art. ID LM24-OP-4011634
Unit 1 mL
Deliverydetails De Minimis UN1888 III E1

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Decane	[124-18-5]	500	mg/L		
	n-Docosane	[629-97-0]	500	mg/L		
	n-Dodecane	[112-40-3]	500	mg/L		
	n-Dotriacontane	[544-85-4]	500	mg/L		
	n-Eicosane	[112-95-8]	500	mg/L		
	n-Heptadecane	[629-78-7]	500	mg/L		
	n-Hexacosane	[630-01-3]	500	mg/L		
	n-Hexadecane	[544-76-3]	500	mg/L		
	n-Hexatriacontane	[630-06-8]	500	mg/L		
	n-Nonane	[111-84-2]	500	mg/L		
	n-Octadecane	[593-45-3]	500	mg/L		
	n-Octane	[111-65-9]	500	mg/L		
	n-Octatriacontane	[7194-85-6]	500	mg/L		
	n-Pentadecane	[629-62-9]	500	mg/L		
	n-Tetradecane	[629-59-4]	500	mg/L		
	n-Tetratriacontane	[14167-59-0]	500	mg/L		
	n-Tridecane	[629-50-5]	500	mg/L		
	n-Undecane	[1120-21-4]	500	mg/L		
	n-Octacosane	[630-02-4]	500	mg/L		
	Phytane	[638-36-8]	500	mg/L		
	Pristane	[1921-70-6]	500	mg/L		
	n-Tetracontane	[4181-95-7]	500	mg/L		
	n-Tetracosane	[646-31-1]	500	mg/L		
	n-Triacontane	[638-68-6]	500	mg/L		