

ASTM Method D3798 p-Xylene Impurity Standard with IS (set of 5x1 mL): 11 components

Art. ID LM24-OP-4011640

Unit 5x1 mL

Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	0,15	%		
	Ethylbenzene	[100-41-4]	0,15	%		
	Isopropylbenzene	[98-82-8]	0,15	%		
	m-Xylene	[108-38-3]	0,15	%		
	n-Octane	[111-65-9]	0,15	%		
	n-Pentane	[109-66-0]	0,15	%		
	n-Undecane	[1120-21-4]	0,5	%		
	o-Xylene	[95-47-6]	0,15	%		
	p-Xylene	[106-42-3]	98,15	%		
	n-Propylbenzene	[103-65-1]	0,15	%		
	Toluene	[108-88-3]	0,15	%		