

## **ASTM Method D5580 Daily Quality Control Standard without IS: 14 components**

Art. ID LM24-OP-4011643

Unit 10 mL

### Description

Validity: 12 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2,4,5-Tetramethylbenzene	[95-93-2]	1	%		
	ene					
	1,2,4-Trimethylbenzene	[95-63-6]	3	%		
	iso-Octane	[540-84-1]	20	%		
	Benzene	[71-43-2]	1	%		
	Ethylbenzene	[100-41-4]	2	%		
	n-Decane	[124-18-5]	10	%		
	n-Heptane	[142-82-5]	20	%		
	n-Hexane	[110-54-3]	12	%		
	n-Octane	[111-65-9]	15	%		
	n-Tridecane	[629-50-5]	1	%		
	Naphthalene	[91-20-3]	1	%		
	o-Xylene	[95-47-6]	2	%		
	p-Xylene	[106-42-3]	3	%		
	Toluene	[108-88-3]	9	%		